

Product Datasheet

Mouse EGF Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit (orb1202792)

| | |
|--------------------------|---------------------------------------------------------------|
| Description | Mouse EGF Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit |
| Reactivity | Mouse |
| Range | 31.2-2000pg/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates or other biological fluids. |
| Uniprot ID | Q9QXT5 |
| Sensitivity | 15.5pg/mL |
| Entrez | 353156 |
| Expiration Date | Please enquire. |