

Product Datasheet

Human Wingless Type MMTV Integration Site Family, Member 5B (WNT5B) ELISA Kit (orb1202844)

Description	Human Wingless Type MMTV Integration Site Family, Member 5B (WNT5B) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H1J7
Sensitivity	0.53ng/mL
Entrez	81029
Expiration Date	Please enquire.