

Product Datasheet

Human Suppressor Of Zeste 12 Homolog (SUZ12) ELISA Kit (orb1202874)

Description	Human Suppressor Of Zeste 12 Homolog (SUZ12) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q15022
Sensitivity	0.113ng/mL
Entrez	23512
Expiration Date	Please enquire.