

Product Datasheet

Human RAB37, Member RAS Oncogene Family (RAB37) ELISA Kit (orb1202916)

Description	Human RAB37, Member RAS Oncogene Family (RAB37) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q96AX2
Sensitivity	0.095ng/mL
Entrez	326624
Expiration Date	Please enquire.