

Product Datasheet

Human Notch2 N-Terminal Like Protein (NOTCH2NL) ELISA Kit (orb1202951)

Description	Human Notch2 N-Terminal Like Protein (NOTCH2NL) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q7Z3S9
Sensitivity	0.117ng/mL
Entrez	388677
Expiration Date	Please enquire.