

Product Datasheet

Human JNK1/MAPK8 Associated Membrane Protein (JKAMP) ELISA Kit (orb1202993)

Description	Human JNK1/MAPK8 Associated Membrane Protein (JKAMP) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9P055
Sensitivity	0.35ng/mL
Entrez	51528
Expiration Date	Please enquire.