

Product Datasheet

Human Interferon Gamma Inducible Protein 16 (IFI16) ELISA Kit (orb1202998)

Description	Human Interferon Gamma Inducible Protein 16 (IFI16) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q16666
Sensitivity	0.049ng/mL
Entrez	3428
Expiration Date	Please enquire.