

Product Datasheet

Human Ghrelin-O-Acyltransferase (GOAT) ELISA Kit (orb1203010)

Description	Human Ghrelin-O-Acyltransferase (GOAT) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q96T53
Sensitivity	0.26ng/mL
Entrez	619373
Expiration Date	Please enquire.