

## Product Datasheet

# Human Fc Fragment Of IgG Binding Protein (FcγBP) ELISA Kit (orb1203017)

<b>Description</b>	Human Fc Fragment Of IgG Binding Protein (FcγBP) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9Y6R7</a>
<b>Sensitivity</b>	0.116ng/mL
<b>Entrez</b>	<a href="#">8857</a>
<b>Expiration Date</b>	Please enquire.