

## Product Datasheet

# Human Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) ELISA Kit (orb1203097)

<b>Description</b>	Human Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	125-8000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q12904</a>
<b>Sensitivity</b>	43pg/mL
<b>Entrez</b>	<a href="#">9255</a>
<b>Expiration Date</b>	Please enquire.