

## Product Datasheet

# Rat VEGF Co Regulated Chemokine 1 (VCC1) ELISA Kit (orb1203113)

<b>Description</b>	Rat VEGF Co Regulated Chemokine 1 (VCC1) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	15.6-1000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">D4A875</a>
<b>Sensitivity</b>	5.7pg/mL
<b>Entrez</b>	<a href="#">308436</a>
<b>Expiration Date</b>	Please enquire.