

Product Datasheet

Rat Sphingosine-1-Phosphate Phosphatase 2 (SGPP2) ELISA Kit (orb1203172)

Description	Rat Sphingosine-1-Phosphate Phosphatase 2 (SGPP2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	F1LZ44
Sensitivity	0.063ng/mL
Entrez	301543
Expiration Date	Please enquire.