

Product Datasheet

Rat Secreted Frizzled Related Protein 4 (SFRP4) ELISA Kit (orb1203177)

Description	Rat Secreted Frizzled Related Protein 4 (SFRP4) ELISA Kit
Reactivity	Rat
Range	0.187-12ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9JLS4
Sensitivity	0.06ng/mL
Entrez	89803
Expiration Date	Please enquire.