

Product Datasheet

Rat Succinate Dehydrogenase Complex Subunit B (SDHB) ELISA Kit (orb1203181)

Description	Rat Succinate Dehydrogenase Complex Subunit B (SDHB) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P21913
Sensitivity	0.041ng/mL
Entrez	298596
Expiration Date	Please enquire.