

Product Datasheet

Rat Amiloride Sensitive Sodium Channel Subunit Alpha (SCNN1a) ELISA Kit (orb1203183)

Description	Rat Amiloride Sensitive Sodium Channel Subunit Alpha (SCNN1a) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P37089
Sensitivity	0.059ng/mL
Entrez	25122
Expiration Date	Please enquire.