

## Product Datasheet

# Rat Protein Phosphatase 2A Activator, Regulatory Subunit 4 (PPP2R4) ELISA Kit (orb1203208)

<b>Description</b>	Rat Protein Phosphatase 2A Activator, Regulatory Subunit 4 (PPP2R4) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>B2RYQ2</b>
<b>Sensitivity</b>	0.057ng/mL
<b>Entrez</b>	<b>362102</b>
<b>Expiration Date</b>	Please enquire.