

Product Datasheet

Rat Neuropeptide Y Receptor Y1 (NPY1R) ELISA Kit (orb1203247)

Description	Rat Neuropeptide Y Receptor Y1 (NPY1R) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P21555
Sensitivity	0.128ng/mL
Entrez	29358
Expiration Date	Please enquire.