

Product Datasheet

Rat Natriuretic Peptide Receptor 2 (NPR2) ELISA Kit (orb1203249)

Description	Rat Natriuretic Peptide Receptor 2 (NPR2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P16067
Sensitivity	0.056ng/mL
Entrez	116564
Expiration Date	Please enquire.