

Product Datasheet

Rat Low Density Lipoprotein Receptor Related Protein 8 (LRP8) ELISA Kit (orb1203281)

Description	Rat Low Density Lipoprotein Receptor Related Protein 8 (LRP8) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	D3ZE75
Sensitivity	13.3pg/mL
Entrez	362558
Expiration Date	Please enquire.