

## Product Datasheet

# Rat Insulin Like Protein 3 (INSL3) ELISA Kit (orb1203300)

<b>Description</b>	Rat Insulin Like Protein 3 (INSL3) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.234-15ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9WUK0</a>
<b>Sensitivity</b>	0.062ng/mL
<b>Entrez</b>	<a href="#">114215</a>
<b>Expiration Date</b>	Please enquire.