

## Product Datasheet

### Canine Clusterin (CLU) ELISA Kit (orb1203418)

<b>Description</b>	Canine Clusterin (CLU) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	1.56-100ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P25473</a>
<b>Sensitivity</b>	0.5ng/mL
<b>Entrez</b>	<a href="#">442971</a>
<b>Expiration Date</b>	Please enquire.