

Product Datasheet

Rat Aldolase C, Fructose Bisphosphate (ALDOC) ELISA Kit (orb1203465)

Description	Rat Aldolase C, Fructose Bisphosphate (ALDOC) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P09117
Sensitivity	0.283ng/mL
Entrez	24191
Expiration Date	Please enquire.