

Product Datasheet

Rat Adrenocorticotrophic Hormone Receptor (ACTHR) ELISA Kit (orb1203481)

Description	Rat Adrenocorticotrophic Hormone Receptor (ACTHR) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	G3V848
Sensitivity	0.234ng/mL
Entrez	282839
Expiration Date	Please enquire.