

Product Datasheet

Bovine Vasoactive Intestinal Peptide (VIP) ELISA Kit (orb1203488)

Description	Bovine Vasoactive Intestinal Peptide (VIP) ELISA Kit
Reactivity	Bovine
Range	7.81-500pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P81401
Sensitivity	3.8pg/mL
Entrez	280956
Expiration Date	Please enquire.