

## Product Datasheet

### Bovine Transthyretin (TTR) ELISA Kit (orb1203497)

<b>Description</b>	Bovine Transthyretin (TTR) ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	7.81-500ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O46375</a>
<b>Sensitivity</b>	2.8ng/mL
<b>Entrez</b>	<a href="#">280948</a>
<b>Expiration Date</b>	Please enquire.