

Product Datasheet

Porcine Surfactant Associated Protein A (SPA) ELISA Kit (orb1203521)

Description	Porcine Surfactant Associated Protein A (SPA) ELISA Kit
Reactivity	Porcine
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, lung lavage fluid or other biological fluids.
Uniprot ID	P49874
Sensitivity	1.16ng/mL
Entrez	397503
Expiration Date	Please enquire.