

Product Datasheet

Porcine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit (orb1203535)

Description	Porcine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
Reactivity	Porcine
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q2EN89
Sensitivity	1.15ng/mL
Entrez	396613
Expiration Date	Please enquire.