

Product Datasheet

Bovine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit (orb1203536)

Description	Bovine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
Reactivity	Bovine
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O97919
Sensitivity	18.8pg/mL
Entrez	327712
Expiration Date	Please enquire.