

## Product Datasheet

### Porcine Osteopontin (OPN) ELISA Kit (orb1203553)

**Description** Porcine Osteopontin (OPN) ELISA Kit

**Reactivity** Porcine

**Range** 62.5-4000pg/mL

**Note** For research use only

**Application notes** sandwich

**Sample Types** serum, plasma, tissue homogenates, cell lysates or other biological fluids.

**Uniprot ID** [P14287](#)

**Sensitivity** 29.8pg/mL

**Entrez** [397087](#)

**Expiration Date** Please enquire.