

Product Datasheet

Porcine Nitric Oxide Synthase 3, Endothelial (NOS3) ELISA Kit (orb1203560)

Description	Porcine Nitric Oxide Synthase 3, Endothelial (NOS3) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q28969
Sensitivity	5.4pg/mL
Entrez	397557
Expiration Date	Please enquire.