

Product Datasheet

Porcine Myeloperoxidase (MPO) ELISA Kit (orb1203563)

Description	Porcine Myeloperoxidase (MPO) ELISA Kit
Reactivity	Porcine
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	K7GRV6
Sensitivity	0.35ng/mL
Entrez	100517120
Expiration Date	Please enquire.