

Product Datasheet

Porcine Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit (orb1203578)

Description	Porcine Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q711P4
Sensitivity	12.7pg/mL
Entrez	396668
Expiration Date	Please enquire.