

Product Datasheet

Bovine Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit (orb1203579)

| | |
|--------------------------|--|
| Description | Bovine Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit |
| Reactivity | Bovine |
| Range | 31.2-2000pg/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | Q17QA1 |
| Sensitivity | 9.8pg/mL |
| Entrez | 414347 |
| Expiration Date | Please enquire. |