

Product Datasheet

Bovine Macrophage Inflammatory Protein 1 Alpha (MIP1a) ELISA Kit (orb1203580)

Description	Bovine Macrophage Inflammatory Protein 1 Alpha (MIP1a) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8SQA6
Sensitivity	7.4pg/mL
Entrez	282170
Expiration Date	Please enquire.