

## Product Datasheet

# Porcine Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit (orb1203581)

<b>Description</b>	Porcine Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	78.1-5000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">C7F842</a>
<b>Sensitivity</b>	32pg/mL
<b>Entrez</b>	<a href="#">100302703</a>
<b>Expiration Date</b>	Please enquire.