

Product Datasheet

Bovine Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit (orb1203582)

Description	Bovine Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q09141
Sensitivity	6.4pg/mL
Entrez	281044
Expiration Date	Please enquire.