

Product Datasheet

Porcine Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit (orb1203583)

Description	Porcine Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit
Reactivity	Porcine
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42831
Sensitivity	0.124ng/mL
Entrez	397422
Expiration Date	Please enquire.