

Product Datasheet

Porcine Hypoxia Inducible Factor 1 Alpha (HIF1a) ELISA Kit (orb1203605)

Description	Porcine Hypoxia Inducible Factor 1 Alpha (HIF1a) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	A0T1H6
Sensitivity	5.8pg/mL
Entrez	396696
Expiration Date	Please enquire.