

## Product Datasheet

### Porcine Cystatin 3 (CST3) ELISA Kit (orb1203640)

<b>Description</b>	Porcine Cystatin 3 (CST3) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, urine or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q0Z8R0</a>
<b>Sensitivity</b>	0.056ng/mL
<b>Entrez</b>	<a href="#">733672</a>
<b>Expiration Date</b>	Please enquire.