

Product Datasheet

Porcine Ciliary Neurotrophic Factor (CNTF) ELISA Kit (orb1203644)

Description	Porcine Ciliary Neurotrophic Factor (CNTF) ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	F1RMF2
Sensitivity	22.4pg/mL
Entrez	100518821
Expiration Date	Please enquire.