

Product Datasheet

Porcine Cholecystokinin (CCK) ELISA Kit (orb1203652)

Description	Porcine Cholecystokinin (CCK) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01356
Sensitivity	4.88pg/mL
Entrez	397468
Expiration Date	Please enquire.