

Product Datasheet

Bovine Beta-Endorphin (bEP) ELISA Kit (orb1203662)

Description	Bovine Beta-Endorphin (bEP) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01190
Sensitivity	4.89pg/mL
Entrez	281416
Expiration Date	Please enquire.