

Product Datasheet

Bovine Aprotinin (AP) ELISA Kit (orb1203668)

Description	Bovine Aprotinin (AP) ELISA Kit
Reactivity	Bovine
Range	1.56-100ng/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P00974
Sensitivity	0.43ng/mL
Entrez	404172
Expiration Date	Please enquire.