

Product Datasheet

Bovine Annexin A5 (ANXA5) ELISA Kit (orb1203669)

Description	Bovine Annexin A5 (ANXA5) ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P81287
Sensitivity	0.16ng/mL
Entrez	281626
Expiration Date	Please enquire.