

Product Datasheet

Porcine Alpha-2-Heremans Schmid Glycoprotein (aHSG) ELISA Kit (orb1203679)

Description	Porcine Alpha-2-Heremans Schmid Glycoprotein (aHSG) ELISA Kit
Reactivity	Porcine
Range	15.6-1000ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P29700
Sensitivity	6.7ng/mL
Expiration Date	Please enquire.