

Product Datasheet

Porcine Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit (orb1203688)

| | |
|--------------------------|---|
| Description | Porcine Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit |
| Reactivity | Porcine |
| Range | 0.39-25ng/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma or other biological fluids. |
| Uniprot ID | F1SN68 |
| Sensitivity | 0.19ng/mL |
| Entrez | 396901 |
| Expiration Date | Please enquire. |