

Product Datasheet

Porcine Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit (orb1203690)

Description	Porcine Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit
Reactivity	Porcine
Range	46.8-3000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P20034
Sensitivity	13.45pg/mL
Entrez	100126843
Expiration Date	Please enquire.