

## Product Datasheet

# Bovine Interleukin 22 (IL22) ELISA Kit (orb1203692)

<b>Description</b>	Bovine Interleukin 22 (IL22) ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	15.6-1000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">A5JS75</a>
<b>Sensitivity</b>	5.4pg/mL
<b>Entrez</b>	<a href="#">507778</a>
<b>Expiration Date</b>	Please enquire.