

## Product Datasheet

# Rat Chromogranin A (CHGA) ELISA Kit (orb1203852)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Rat Chromogranin A<br>(CHGA) ELISA Kit       |
| <b>Reactivity</b>        | Rat  |
| <b>Range</b>             | 31.2-2000pg/mL                               |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | sandwich                                     |
| <b>Sample Types</b>      | serum, plasma or other<br>biological fluids. |
| <b>Uniprot ID</b>        | <a href="#">Q5PPG5</a>                       |
| <b>Sensitivity</b>       | 12.3pg/mL                                    |
| <b>Entrez</b>            | <a href="#">24258</a>                        |
| <b>Expiration Date</b>   | Please enquire.                              |