

Product Datasheet

Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit (orb1203907)

Description	Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q8R4B8
Sensitivity	0.055ng/mL
Entrez	216799
Expiration Date	Please enquire.